

Anti-PMS2 antibody [IHC422] (STJA0042248)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-p
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide
Isotype	IgGk (Rab)
Storage Instruction	Store at 2-8°C for up to 1 year. Do not freeze.

TARGET INFORMATION

Gene ID	5395
Gene Symbol	PMS2
UniProt ID	PMS2_HUMAN
Specificity	Tested and validated on colon cancer tissue.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---