

Anti-beta Tubulin antibody [2E8-2] {FITC} (STJA0042294)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	FITC
Purification	Protein A/G Purification
Dilution Range	FCM 0.2 μ g per 10^6 cells in a 100 μ L suspension
Formulation	PBS with 0.02% sodium azide and 0.1% gelatin, pH7.4
Isotype	IgG
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	203068
Gene Symbol	TUBB
UniProt ID	TBB5_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---