

## Anti-alpha Tubulin antibody [2H2] {EF647} (STJA0042298)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	EF647
Purification	Protein A/G Purification
Dilution Range	FCM 0.2 $\mu$ g per $10^6$ cells in a 100 $\mu$ L suspension
Formulation	PBS with 0.02% sodium azide, 1% protective protein and 50% glycerol, pH7.4
Isotype	IgG2a
Storage Instruction	Store at -20 °C Valid for 12 months. Avoid freeze/thaw cycles. Protected from prolonged exposure to light.

### TARGET INFORMATION

Gene ID	<a href="#">10376</a>
Gene Symbol	<a href="#">TUBA1B</a>
UniProt ID	<a href="#">TBA1B_HUMAN</a>

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---