

## Anti-GAPDH antibody [7B8] {HRP} (STJA0042301)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human/Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	HRP
Purification	Protein A/G Purification
Dilution Range	WB 1:2500-1:5000
Formulation	PBS with 0.1% procline 300, 1% protective protein and 50% glycerol, pH7.4
Isotype	IgG1
Storage Instruction	Store at-20 °C Valid for 12 months. Avoid freeze/thaw cycles. Protected from prolonged exposure to light.

### TARGET INFORMATION

Gene ID	<a href="#">2597</a>
Gene Symbol	<a href="#">GAPDH</a>
UniProt ID	<a href="#">G3P_HUMAN</a>

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---