

Anti-GAPDH antibody [7B8] {EF594} (STJA0042303)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	EF594
Purification	Protein A/G Purification
Dilution Range	IF 1:100-1:200
Formulation	PBS with 0.02% sodium azide, 1% protective protein and 50% glycerol, pH7.4
Isotype	IgG1
Storage Instruction	Store at-20 °C Valid for 12 months. Avoid freeze/thaw cycles. Protected from prolonged exposure to light.

TARGET INFORMATION

Gene ID	2597
Gene Symbol	GAPDH
UniProt ID	G3P_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---