

Anti-Vinculin antibody [4D4] {EF647} (STJA0042307)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	EF647
Purification	Protein A/G Purification
Dilution Range	FCM 0.2 µg per 10 ⁶ cells in a 100 µL suspension
Formulation	PBS with 0.02% sodium azide, 1% protective protein and 50% glycerol, pH7.4
Isotype	IgG1
Storage Instruction	Store at-20 °C Valid for 12 months. Avoid freeze/thaw cycles. Protected from prolonged exposure to light.

TARGET INFORMATION

Gene ID	7414
Gene Symbol	VCL
UniProt ID	VINC_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---