

Anti-Lamin A/C antibody [4F1] {EF488} (STJA0042310)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IF |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 1 mg/mL |
| Conjugation | EF488 |
| Purification | Protein A/G Purification |
| Dilution Range | IF 1:100-1:200 |
| Formulation | PBS with 0.02% sodium azide, 1% protective protein and 50% glycerol, pH7.4 |
| Isotype | IgG1 |
| Storage Instruction | Store at -20 °C Valid for 12 months. Avoid freeze/thaw cycles. Protected from prolonged exposure to light. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 4000 |
| Gene Symbol | LMNA |
| UniProt ID | LMNA_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|