

## Anti-Recombinant-Flag-Tag antibody [3H9] {EF594} (STJA0042314)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF
Host / Source	Mouse
Reactivity	All

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	EF594
Purification	Protein A/G Purification
Dilution Range	IF 1:100-1:200
Formulation	PBS with 0.02% sodium azide, 1% protective protein and 50% glycerol, pH7.4
Isotype	IgG2a
Storage Instruction	Store at -20 °C Valid for 12 months. Avoid freeze/thaw cycles. Protected from prolonged exposure to light.

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---