

Anti-Recombinant-Myc-Tag antibody [6B10] {FITC} (STJA0042315)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	All

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	FITC
Purification	Protein A/G Purification
Dilution Range	FCM 0.2 μ g per 10^6 cells in a 100 μ L suspension
Formulation	PBS with 0.02% sodium azide and 0.1% gelatin, pH7.4
Isotype	IgG1
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---