

Anti-Recombinant-HA-Tag antibody [4F5] {EF647} (STJA0042321)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	All

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	EF647
Purification	Protein A/G Purification
Dilution Range	FCM 0.2 μ g per 10^6 cells in a 100 μ L suspension
Formulation	PBS with 0.02% sodium azide, 1% protective protein and 50% glycerol, pH7.4
Isotype	IgG2a
Storage Instruction	Store at -20 °C Valid for 12 months. Avoid freeze/thaw cycles. Protected from prolonged exposure to light.

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---