

Anti-CD122 antibody [TM-Beta 1] {PE/Cyanine 5} (STJA0042464)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine 5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16185
Gene Symbol	Il2rb
UniProt ID	IL2RB_MOUSE

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---