

## Anti-CD8 antibody [UCHT-4] {PerCP/Cyanine 5.5} (STJA0042523)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 5 Mu L/Test  |
| Conjugation         | PerCP/Cyanine 5.5  |
| Formulation         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.      |
| Isotype             | IgG1k  |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">925</a>        |
| Gene Symbol | <a href="#">CD8A</a>       |
| UniProt ID  | <a href="#">CD8A_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|