

Anti-MHC I antibody [SF1.1.10] {PE} (STJA0042557)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Mouse |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|---------------------|--------------------------------------------------------------------------------------------------|
| Clonality | Monoclonal |
| Concentration | 0.2 mg/mL |
| Conjugation | PE |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG2ak |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 14972 |
| Gene Symbol | H2-K1 |
| UniProt ID | HA1B_MOUSE |

ADDITIONAL INFORMATION

| | |
|------|-------------------------------------------------------------------------------------------------------------------|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------|