

Anti-GFP antibody [7C4] {EF488} (STJA0042570)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF
Host / Source	Mouse
Reactivity	All

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	EF488
Purification	Protein A/G Purification
Dilution Range	IF 1:500-1:2000
Formulation	PBS with 0.02% sodium azide, 1% protective protein and 50% glycerol, pH 7.4
Isotype	IgG2a
Storage Instruction	Store at -20 °C Valid for 12 months. Avoid freeze/thaw cycles. Protected from prolonged exposure to light.

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---