

## Anti-GFP Tag antibody [7C4] {EF594} (STJA0042571)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IF                 |
| Host / Source | Mouse              |
| Reactivity    | All                |

### PRODUCT PROPERTIES

|                     |                                                                                                            |
|---------------------|------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                 |
| Concentration       | 1 mg/mL                                                                                                    |
| Conjugation         | EF594                                                                                                      |
| Purification        | Protein A/G Purification                                                                                   |
| Dilution Range      | IF 1:500-1:2000                                                                                            |
| Formulation         | PBS with 0.02% sodium azide, 1% protective protein and 50% glycerol, pH7.4                                 |
| Isotype             | IgG2a                                                                                                      |
| Storage Instruction | Store at -20 °C Valid for 12 months. Avoid freeze/thaw cycles. Protected from prolonged exposure to light. |

### ADDITIONAL INFORMATION

|      |                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------|