

Anti-PCNA antibody [9C6] {HRP} (STJA0042572)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	HRP
Purification	Protein A/G Purification
Dilution Range	WB 1:2500-1:5000
Formulation	PBS with 0.1% procline 300, 1% protective protein and 50% glycerol, pH7.4
Isotype	IgG2a
Storage Instruction	Store at-20 °C Valid for 12 months. Avoid freeze/thaw cycles. Protected from prolonged exposure to light.

TARGET INFORMATION

Gene ID	5111
Gene Symbol	PCNA
UniProt ID	PCNA_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---