

Anti-CD8b.2 antibody [53-5.8] {FITC} (STJA0042579)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 5 Mu L/Test |
| Conjugation | FITC |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG1k |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 12526 |
| Gene Symbol | Cd8b |
| UniProt ID | CD8B_MOUSE |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|