

Anti-HRP-AVI-Tag antibody [1A3] {HRP} (STJA0042701)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	All

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	HRP
Purification	Antigen Affinity Purification
Dilution Range	WB 1:4000-20000 IP 3ug/sample IF 1:100-200
Formulation	PBS with 0.1% procline 300, 1% protective protein and 50% glycerol, pH7.4
Isotype	IgG1
Storage Instruction	Store at -20°C for the long term, and avoid freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---