

Anti-CD11c antibody [N418] {EF421} (STJA0042764)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	FC
	Host / Source	Armenian Hamster
	Reactivity	Mouse
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Concentration	5 Mu L/Test
	Conjugation	EF421
	Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
	Isotype	IgG2
	Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
TARGET INFORMATION		
	Gene ID	16411
	Gene Symbol	Itgax
	UniProt ID	ITAX_MOUSE
ADDITIONAL INFORMATION		
	Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.