

Anti-IL-4 antibody [MP4-25D2] {EF700} (STJA0043004)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	EF700
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3565
Gene Symbol	IL4
UniProt ID	IL4_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---