

**Anti-CD16 antibody [3G8] {EF540} (STJA0043036)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human/Rhesus/Cynomolgus

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	EF540
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

**TARGET INFORMATION**

Gene ID	<a href="#">2215</a> <a href="#">2214</a>
Gene Symbol	<a href="#">FCGR3B</a> <a href="#">FCGR3A</a>
UniProt ID	<a href="#">FCG3B_HUMAN</a> <a href="#">FCG3A_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---