

## Human ACE (Angiotensin-Converting Enzyme) Antibody Pair (STJAP00008)

### GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Goat/Goat
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Storage Instruction	Store at -20°C for the long term, and avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">1636</a>
Gene Symbol	<a href="#">ACE</a>
UniProt ID	<a href="#">ACE_HUMAN</a>

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---