

Rat IgG (Ig Gamma-2c Chain C Region) Antibody Pair (STJAP00090)

GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Rabbit/Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

UniProt ID	IGG2C_RAT
------------	---------------------------

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---