

Rat KIM-1 (Hepatitis A Virus Cellular Receptor 1 Homolog) Antibody Pair (STJAP00096)

GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Mouse/Goat
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal/Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Antigen Affinity.
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Storage Instruction	Store at -20°C for the long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	286934
Gene Symbol	Havcr1
UniProt ID	HAVR1_RAT

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---