

Human HSP-70/HSPA9 (Stress-70 Protein-Mitochondrial) Antibody Pair (STJAP00175)

GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Rabbit/Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A & Antigen Affinity. Detection Purification: Protein A.
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3313
Gene Symbol	HSPA9
UniProt ID	GRP75_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---