

Human PGF (Placenta Growth Factor) Antibody Pair (STJAP00205)

GENERAL INFORMATION

| | |
|---------------|----------------|
| Product Type | Antibody Pairs |
| Applications | ELISA |
| Host / Source | Mouse/Goat |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal/Polyclonal |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated/Biotin |
| Purification | Capture Purification: Protein A or G. Detection Purification: Antigen Affinity. |
| Formulation | Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. |
| Storage Instruction | Store at -20°C for the long term, and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 5228 |
| Gene Symbol | PGF |
| UniProt ID | PLGF_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|