

## Human PDGF-AB (PDGF-AB) Antibody Pair (STJAP00206)

### GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Mouse/Goat
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal/Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Antigen Affinity.
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Storage Instruction	Store at -20°C for the long term, and avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">5154</a> <a href="#">5155</a>
Gene Symbol	<a href="#">PDGFA</a> <a href="#">PDGFB</a>
UniProt ID	<a href="#">PDGFA_HUMAN</a> <a href="#">PDGFB_HUMAN</a>

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---