

Human FG (Fibrinogen Alpha Chain Cleaved Into-Fibrinopeptide A-Fibrinogen Alpha Chain) Antibody Pair (STJAP00236)

GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Sheep/Sheep
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Storage Instruction	Store at -20°C for the long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2243
Gene Symbol	FGA
UniProt ID	FIBA_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---