

Mouse CTLA4 (Cytotoxic T-Lymphocyte Protein 4) Antibody Pair (STJAP00306)

GENERAL INFORMATION

| | |
|---------------|----------------|
| Product Type | Antibody Pairs |
| Applications | ELISA |
| Host / Source | Goat/Goat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Polyclonal/Polyclonal |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated/Biotin |
| Purification | Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity. |
| Formulation | Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. |
| Storage Instruction | Store at -20°C for the long term, and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 12477 |
| Gene Symbol | Ctla4 |
| UniProt ID | CTLA4_MOUSE |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|