

## Mouse IGF-2 (Insulin-Like Growth Factor 2) Antibody Pair (STJAP00314)

### GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Goat/Rat
Reactivity	Mouse

### PRODUCT PROPERTIES

Clonality	Polyclonal/Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Protein A or G.
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Storage Instruction	Store at -20°C for the long term, and avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">16002</a>
Gene Symbol	<a href="#">Igf2</a>
UniProt ID	<a href="#">IGF2_MOUSE</a>

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---