

Rat Beta TG/PBP/CXCL7/NAP2 (Beta TG/PBP/CXCL7/NAP2) Antibody Pair (STJAP00357)**GENERAL INFORMATION**

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Goat/Goat
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal/Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity.
Formulation	Capture Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4. Detection Buffer: PBS with 0.04% Proclin 300, 50% glycerol, pH 7.4.
Storage Instruction	Store at -20°C for the long term, and avoid freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---