

Human TACE/ADAM17 (Disintegrin And Metalloproteinase Domain-Containing Protein 17) Antibody Pair (STJAP00444)

GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Mouse/Mouse
Reactivity	Human

PRODUCT PROPERTIES

Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Protein A or G.
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Storage Instruction	Store at -20°C for the long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6868
Gene Symbol	ADAM17
UniProt ID	ADA17_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---