

Human RNASE3/ECP (Eosinophil Cationic Protein) Antibody Pair (STJAP00457)

GENERAL INFORMATION

| | |
|---------------|----------------|
| Product Type | Antibody Pairs |
| Applications | ELISA |
| Host / Source | Rabbit/Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated/Biotin |
| Purification | Capture Purification: Antigen Affinity. Detection Purification: Antigen Affinity. |
| Formulation | Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5. |
| Storage Instruction | Store at -20°C for the long term, and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------|---------------------------|
| Gene ID | 6037 |
| Gene Symbol | RNASE3 |
| UniProt ID | ECP_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|