

## Porcine IL-1 Alpha (IL-1 Alpha) Antibody Pair (STJAP00587)

### GENERAL INFORMATION

|               |                |
|---------------|----------------|
| Product Type  | Antibody Pairs |
| Applications  | ELISA          |
| Host / Source | Goat/Goat      |
| Reactivity    | Porcine        |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Concentration       | 0.5 mg/mL   |
| Conjugation         | Unconjugated/Biotin   |
| Purification        | Capture Purification: Antigen Affinity.<br>Detection Purification: Antigen Affinity.  |
| Formulation         | Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5.<br>Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5. |
| Storage Instruction | Store at -20°C for the long term, and avoid freeze-thaw cycles.   |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|