

## Rat RBP4 (Retinol-Binding Protein 4) Antibody Pair (STJAP00609)

### GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Rabbit/Rabbit
Reactivity	Rat

### PRODUCT PROPERTIES

Concentration	0.5 mg/mL
Conjugation	Unconjugated/Biotin
Purification	Capture Purification: Protein A or G. Detection Purification: Antigen Affinity.
Formulation	Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5. Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5.
Storage Instruction	Store at -20°C for the long term, and avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">25703</a>
Gene Symbol	<a href="#">Rbp4</a>
UniProt ID	<a href="#">RET4_RAT</a>

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---