

## Monkey RETN (RETN) Antibody Pair (STJAP00632)

### GENERAL INFORMATION

|               |                |
|---------------|----------------|
| Product Type  | Antibody Pairs |
| Applications  | ELISA          |
| Host / Source | Mouse/Mouse    |
| Reactivity    | Monkey         |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                 |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Concentration       | 0.5 mg/mL                                                                                                                                                       |
| Conjugation         | Unconjugated/Biotin                                                                                                                                             |
| Purification        | Capture Purification: Protein A or G.<br>Detection Purification: Protein A or G.                                                                                |
| Formulation         | Capture Buffer: PBS with 0.04% Proclin 300; 50% glycerol; pH 7.5.<br>Detection Buffer: PBS with 0.04% Proclin 300; 1% protective protein; 50% glycerol; pH 7.5. |
| Storage Instruction | Store at -20°C for the long term, and avoid freeze-thaw cycles.                                                                                                 |

### ADDITIONAL INFORMATION

|      |                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------|