

Human TNF alpha (TNF alpha) Antibody Pair (STJAP00661)

GENERAL INFORMATION

| | |
|---------------|----------------|
| Product Type | Antibody Pairs |
| Applications | ELISA |
| Host / Source | Mouse/Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal/Monoclonal |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | Affinity Purified |
| Dilution Range | ELISA 1:1000 |
| Formulation | Liquid in PBS containing 0.09% sodium azide as a preservative./Liquid in PBS containing 0.01% thimerosal as a preservative. |
| Isotype | IgG1/IgG2a |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 7124 |
| Gene Symbol | TNE |
| UniProt ID | TNFA_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|