

Human S100B (S100B) Antibody Pair (STJAP00668)

GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Mouse/Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	ELISA 1:1000
Formulation	Liquid in PBS containing 0.09% sodium azide as a preservative./Liquid in PBS containing 0.01% thimerosal as a preservative.
Isotype	IgG1/IgG2a
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	6285
Gene Symbol	S100B
UniProt ID	S100B_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---