

## Human IL-1 beta (IL-1 beta) Antibody Pair (STJAP00671)

### GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Mouse/Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	ELISA 1:1000
Formulation	Liquid in PBS containing 0.09% sodium azide as a preservative./Liquid in PBS containing 0.01% thimerosal as a preservative.
Isotype	IgG2a/IgG3a
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">3553</a>
Gene Symbol	<a href="#">IL1B</a>
UniProt ID	<a href="#">IL1B_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---