

## Human MCP-1 (MCP-1) Antibody Pair (STJAP00675)

### GENERAL INFORMATION

|               |                |
|---------------|----------------|
| Product Type  | Antibody Pairs |
| Applications  | ELISA          |
| Host / Source | Mouse/Mouse    |
| Reactivity    | Human          |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal/Monoclonal   |
| Concentration       | 1 mg/mL   |
| Conjugation         | Unconjugated  |
| Purification        | Affinity Purified   |
| Dilution Range      | ELISA 1:1000  |
| Formulation         | Liquid in PBS containing 0.09% sodium azide as a preservative./Liquid in PBS containing 0.01% thimerosal as a preservative. |
| Isotype             | IgG2b   |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                    |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">6347</a>       |
| Gene Symbol | <a href="#">CCL2</a>       |
| UniProt ID  | <a href="#">CCL2_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|