

Rat IL-10 (IL-10) Antibody Pair (STJAP00751)

GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Rabbit/Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	Determined by end user
Formulation	0.01M PBS, pH 7.2
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	25325
Gene Symbol	Il10
UniProt ID	IL10_RAT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---