

Mouse MMP10 (MMP10) Antibody Pair (STJAP00783)

GENERAL INFORMATION

| | |
|---------------|----------------|
| Product Type | Antibody Pairs |
| Applications | ELISA |
| Host / Source | Rabbit/Rabbit |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal/Monoclonal |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | Affinity Purified |
| Dilution Range | Determined by end user |
| Formulation | 0.01M PBS, pH 7.2 |
| Isotype | IgG |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 17384 |
| Gene Symbol | Mmp10 |
| UniProt ID | MMP10_MOUSE |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|