

Human Amyloid Precursor Protein (Amyloid Precursor Protein) Antibody Pair (STJAP00844)

GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Rabbit/Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal/Monoclonal
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	Determined by end user
Formulation	0.01M PBS, pH 7.2
Isotype	IgG
Storage Instruction	Store at 2-8°C. Do not freeze.

TARGET INFORMATION

Gene ID	351
Gene Symbol	APP
UniProt ID	A4_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---