

## Human ErbB2/HER2 (ErbB2/HER2) Antibody Pair (STJAP00849)

### GENERAL INFORMATION

|               |                |
|---------------|----------------|
| Product Type  | Antibody Pairs |
| Applications  | ELISA          |
| Host / Source | Rabbit/Rabbit  |
| Reactivity    | Human          |

### PRODUCT PROPERTIES

|                     |                                |
|---------------------|--------------------------------|
| Clonality           | Monoclonal/Monoclonal          |
| Conjugation         | Unconjugated                   |
| Purification        | Affinity Purified              |
| Dilution Range      | Determined by end user         |
| Formulation         | 0.01M PBS, pH 7.2              |
| Isotype             | IgG                            |
| Storage Instruction | Store at 2-8°C. Do not freeze. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">2064</a>        |
| Gene Symbol | <a href="#">ERBB2</a>       |
| UniProt ID  | <a href="#">ERBB2_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|