

## Human IDH1 (IDH1) Antibody Pair (STJAP00878)

### GENERAL INFORMATION

|               |                |
|---------------|----------------|
| Product Type  | Antibody Pairs |
| Applications  | ELISA          |
| Host / Source | Rabbit/Rabbit  |
| Reactivity    | Human          |

### PRODUCT PROPERTIES

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal/Monoclonal                                                                    |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purification                                                                    |
| Dilution Range      | Determined by end user                                                                   |
| Formulation         | 0.01M PBS, pH 7.2                                                                        |
| Isotype             | IgG                                                                                      |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">3417</a>       |
| Gene Symbol | <a href="#">IDH1</a>       |
| UniProt ID  | <a href="#">IDHC_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|