

## Human TSP-4 (TSP-4) Antibody Pair (STJAP00939)

### GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Goat/Goat
Reactivity	Human

### PRODUCT PROPERTIES

Concentration	0.5 mg/mL/0.5-4ug/mL
Conjugation	Unconjugated/Biotin
Dilution Range	1:1000-1:10000
Storage Instruction	Store at-20°C for the long term, and avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">7060</a>
Gene Symbol	<a href="#">THBS4</a>
UniProt ID	<a href="#">TSP4_HUMAN</a>

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---