

## Mouse CD117 (CD117) Antibody Pair (STJAP00954)

### GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Rabbit/Rabbit
Reactivity	Mouse

### PRODUCT PROPERTIES

Concentration	0.5 mg/mL/0.5-4ug/mL
Conjugation	Unconjugated/Biotin
Dilution Range	1:1000-1:10000
Storage Instruction	Store at-20°C for the long term, and avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">16590</a>
Gene Symbol	<a href="#">Kit</a>
UniProt ID	<a href="#">KIT_MOUSE</a>

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---