

Human IL-15 (IL-15) Antibody Pair (STJAP00960)

GENERAL INFORMATION

| | |
|---------------|----------------|
| Product Type | Antibody Pairs |
| Applications | ELISA |
| Host / Source | Mouse/Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Concentration | 0.5 mg/mL/0.5-4ug/mL |
| Conjugation | Unconjugated/Biotin |
| Dilution Range | 1:1000-1:10000 |
| Storage Instruction | Store at-20°C for the long term, and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 3600 |
| Gene Symbol | IL15 |
| UniProt ID | IL15_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS. |
|------|---|