

Mouse NOV (NOV) Antibody Pair (STJAP00961)

GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Rabbit/Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Concentration	0.5 mg/mL/0.5-4ug/mL
Conjugation	Unconjugated/Biotin
Dilution Range	1:1000-1:10000
Storage Instruction	Store at-20°C for the long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	18133
Gene Symbol	Ccn3
UniProt ID	CCN3_MOUSE

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---