

Human FGF1 (FGF1) Antibody Pair (STJAP00963)

GENERAL INFORMATION

Product Type	Antibody Pairs
Applications	ELISA
Host / Source	Goat/Goat
Reactivity	Human

PRODUCT PROPERTIES

Concentration	0.5 mg/mL/0.5-4ug/mL
Conjugation	Unconjugated/Biotin
Dilution Range	1:1000-1:10000
Storage Instruction	Store at-20°C for the long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2246
Gene Symbol	FGF1
UniProt ID	FGF1_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---